



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

Trident 10X PBST

Cat No. GTX30977

Application WB, ELISA

Package
100 ml

PRODUCT

Product Description Phosphate-Buffered Saline with Tween Solution (10X)

APPLICATION

Application Note

Dilute this 10X buffer 10 times (e.g. 90 ml of deionized or distilled water + 10 ml of this buffer), mix well, and store 1 X buffer solution at room temperature or 2-8°C for several weeks. Wash membrane or ELISA plates for 10 minutes each for three times. For 10x10 cm membrane use at least 50 ml buffer with shaking gently. This buffer is used in between steps in immunoassays protocol. For diluting blocking solution, primary antibody or HRP-secondary antibody, please refer to WB or ELISA protocol.

PROPERTIES

Form Liquid

Storage Room temperature. If stored at 2-8°C, some sodium chloride will be precipitated, if this occurs, bring the buffer to room temperature and mix well till all sodium chloride is dissolved.

Note For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).