



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Signal+ for Immunostaining

Cat No. GTX49998

Application	ICC/IF, IHC
-------------	-------------

Package	1 kit
---------	-------

PRODUCT

Product Description	Signal+ Antibody Enhancer Solution for Immunostaining (ICC/IF and IHC)
---------------------	--

APPLICATION

Application Note

Signal+ for Immunostaining Solutions A and B exhibit different enhancing effects, depending on antigens and antibodies. These solutions can be used independently. However, both solutions should be explored in order to select the more suitable one for your particular assay system. Solutions are available for purchase individually upon request.


PROPERTIES

Form	Liquid
------	--------

Storage	Store at 4°C.
---------	---------------

Note	For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
------	--

DATA IMAGES




GTX49998 ICC/IF Image


Stathmin antibody diluted in Signal+ stains with an intensity 2-5x higher than antibody diluted in blocking buffer



For full product information, images and publications, please visit our [website](#).


GTx49998 ICC/IF Image

XRCC5 antibody diluted in Signal+ stains with an intensity 2-5x higher than antibody diluted in blocking buffer.


GTx49998 IHC Image

Brain tissue section was stained with secondary antibody conjugated with AF555 (red) and counter stained with Hoechst dye (blue).



For full product information, images and publications, please visit our [website](#).