



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BCIP/NBT chromogen/substrate

Cat No. GTX73326

| | |
|-------------|-----|
| Application | IHC |
|-------------|-----|

| | |
|---------|-------|
| Package | 15 ml |
|---------|-------|

PRODUCT

| | |
|---------------------|---|
| Product Description | Alkaline Phosphatase Substrate-Chromogen (BCIP/NBT) |
|---------------------|---|

APPLICATION

Application Note

1. Once tissue sections have been incubated with alkaline phosphatase, wash them with buffer thoroughly. 2. Wipe the glass to remove excess of buffer and add enough drops of the BCIP/NBT solution to cover the tissue sections. 3. Incubate for 10-30 minutes at room temperature. For the best results, look under the microscope for the signal development. Once desired signal to noise ratio is achieved, stop the reaction by washing the slides in the wash buffer.

PROPERTIES

| | |
|---------|---|
| Storage | Keep as concentrated solution, aliquot and store at 4°C. Do not freeze. |
|---------|---|

| | |
|------|--|
| Note | For laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. |
|------|--|



For full product information, images and publications, please visit our [website](#).