



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

Competent *E. coli* Top10

Cat. No.: 5-1600-020

Version: 9.0
Revision Date: 08.11.2021

Description	Competent <i>E. coli</i> TOP10 cells are ready for heat shock transformation with vector DNA and its subsequent propagation for cloning and transfection purposes.
Genotype	F ⁻ <i>mcrA</i> Δ(<i>mrr-hsdRMS-mcrBC</i>) Φ80 <i>lacZ</i> Δ <i>M15</i> Δ <i>lacX74</i> <i>recA1</i> <i>ara</i> Δ139 Δ(<i>ara-leu</i>)7697 <i>galU</i> <i>galK</i> <i>rpsL</i> (Str ^R) <i>endA1</i> <i>nupG</i>
Transformation efficiency	>1 x 10 ⁷ cfu/μg supercoiled DNA
Form	20 separate one-shot reactions
Amount	100 μl per reaction
Biosafety level	S1 (Germany), BSL-2 (USA)
Stability	12 months after shipping
Storage	Store in cryo storage system at -90 °C to -60 °C
Shipping	Dry ice
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	Thaw a vial of competent <i>E. coli</i> TOP10 on ice. Add up to 10 μl DNA (e.g., from a StarGate® ligation reaction) to the thawed competent <i>E. coli</i> TOP10 cells. Mix gently (do not vortex) and then incubate for 30 min on ice. Repeat gentle mixing and incubate for 5 min at 37 °C. Mix gently and then incubate for 2-5 min on ice. Add 900 μl LB medium and shake for 45 min at 37 °C. This incubation step is necessary especially when using kanamycin to express resistance genes prior to plating on plates for selection. Plate 100 μl on LB agar containing antibiotic (if required) and 50 mg/L X-gal (optional). Centrifuge the residual 900 μl cell mixture for 30 sec in a microfuge, resuspend the cells with 100 μl LB medium and plate the whole amount as above. Incubate plates over night at 37 °C. Pick single colonies for further analyses (plasmid isolation, PCR...).

For research use only

Trademark information

The owners of trademarks marked by "®" or "TM" are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.