



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Anti-Mouse MHC Class II In Vivo Antibody – Low Endotoxin (M5/114.15.2) [ICH1176]

Description

Bulk anti-MHC Class II In Vivo Antibody – Low Endotoxin (M5/114.15.2)

Product Benefits:

ichorbio's anti-MHC Class II In Vivo Antibody – Low Endotoxin (M5/114.15.2) is manufactured in a cGMP compliant, ISO Quality Standard 9001:2015 facility. ichorbio's low endotoxin antibodies have half the endotoxin of comparable antibodies from our [competitors](#) at less than 1.0 EU/mg. If ichorbio's low endotoxin antibodies are not low enough we also offer ultra low endotoxin antibodies which have even less endotoxin (<0.75EU/mg) at an even higher purity (98% versus 95%). ichorbio: the best antibodies for *in vivo* research.

Target:

MHC Class II

Other Names:

I-A/I-E MHC class II

Clone:

M5/114.15.2

Size:

ichorbio's M5/114.15.2 *in vivo* antibody is available in the following bulk sizes: 1mg, 5mg, 25mg, 50mg and 100mg. ichorbio regularly manufactures multi-gram amounts of our anti-MHC Class II In Vivo Antibody – Low Endotoxin (M5/114.15.2) – please [contact us](#) for bulk pricing.

Isotype:

Rat IgG2b

Host:

Rat

Species Reactivity:

Mouse

Specificity:

Clone M5/114.15.2 recognizes an epitope on the mouse MHC class II molecules I-Ab, I-Ad, I-Aq, I-Ed, and I-Ek. Clone M5/114.15.2 does not react with I-Af, I-Ak, or I-As.

Purification Method:

This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Antigen Distribution:

I-Ab/d/q and I-Ed/k are expressed on antigen-presenting cells, including macrophages, monocytes, DCs, and B cells, and activated T cells from mice of the I-Ab/d/q and I-Ed/k haplotypes.

Background:

Anti-MHC Class II In Vivo Antibody – Low Endotoxin (M5/114.15.2), recognizes the major histocompatibility complex (MHC) class II molecules I-Ab, I-Ad, I-Aq, I-Ed, and I-Ek. MHC class II is constitutively expressed on professional antigen-presenting cells (APCs), including macrophages/monocytes, dendritic cells (DCs), and B cells, and is induced on T cells upon activation. MHC class II consists of two transmembrane proteins, a 35 kDa ? (heavy) chain and 29 kDa ? (light) chain. The N-terminal ?1 and ?1 domains form the antigen-binding groove, which binds 13-25 aa peptides derived from exogenous antigens. On APCs, MHC class II plays a critical role in the adaptive immune response by presenting phagocytosed antigens to helper CD4 T cells. The T cell receptor (TCR)/CD3 complex of CD4 T cells interacts with peptide-MHC class II, which induces CD4 T cell activation leading to the coordination and regulation of other effector cells. CD4 molecules also bind to MHC class II, which helps augment TCR signaling. It has also been demonstrated that MHC class II express on activated T cells are capable of antigen presentation and can transduce signals into T cells, enhancing T cell proliferation and activity.

Immunogen:

C57BL/6 Spleen Cells.

Concentration:

>2.0 mg/ml

Formulation:

This M5/114.15.2 *in vivo* antibody is aseptically packaged and formulated in 0.01 M phosphate buffered saline (PBS) pH 7.2 – 7.4, 150 mM NaCl with no carrier protein, potassium or preservatives added. BSA and azide free

Purity:

>95% by SDS-PAGE and HPLC

>98% by SDS-PAGE and HPLC

Endotoxin:

? 1.0 EU/mg as determined by the LAL method

? 0.75 EU/mg as determined by the LAL method

Aggregation:

Aggregation level ? 5%

Aggregation level ? 1%

IMPACT Pathogen Test:

We use the IMPACT test generated by IDEXX Laboratories to guarantee our Ultra Low Endotoxin antibodies are pathogen free. Our rat antibodies are tested for: Mycoplasma spp. Mycoplasma pulmonis Sendai virus Pneumonia virus of mice Rat Minute virus Rat parvovirus Lymphocytic choriomeningitis virus Rat cytomegalovirus Rat coronavirus Sialodacryoadenitis virus Seoul virus Mouse adenovirus Reovirus 3 Rat teilovirus

Storage:

This M5/114.15.2 *in vivo* antibody is stable for at least one week when stored sterile at 2-8°C. For long term storage aseptically aliquot in working volumes without diluting and store at –80°C. Avoid Repeated Freeze Thaw Cycles.

Applications:

IHC (Frozen), Flow Cytometry, Immunoprecipitation, Blocking, CODEX®

Application Notes:

Each investigator should determine their own optimal working dilution for specific applications.

Use:

Products are for research use only. Not for use in diagnostic or therapeutic procedures.