



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



## Cusatuzumab Biosimilar – Research Grade [ICH5037]

### Description

## Cusatuzumab Biosimilar – Research Grade

#### Product Benefits:

ichorbio's Cusatuzumab Biosimilar for Research Use Only (RUO) and not for therapeutic use.

#### Target:

CD70

#### Size:

ichorbio's Cusatuzumab Biosimilar is available in the following sizes:  
1mg, 5mg, 10mg, 20mg, 50mg, 100mg

#### Isotype:

Human IgG1-Lambda2

#### Host:

CHO cells

#### Species Reactivity:

Human

#### Specificity:

Detects human CD70. This non-therapeutic antibody uses the same sequence as the therapeutic antibody Cusatuzumab.

---

Purification Method:

This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

Concentration:

5mg/ml

Formulation:

Sterile, preservative-free, solution in PBS. BSA and Azide free.

Purity:

>95% by SDS-PAGE

Endotoxin:

? 1.0 EU/mg as determined by the LAL method

Storage:

This antibody is stable for at least 4 weeks when stored at 2-8°C. For long term storage, aliquot in working volumes without diluting and store at – 20°C or -80°C. Avoid repeated freeze thaw cycles.

Applications:

Functional Assays

Application Notes:

This biosimilar is strictly for research use only (RUO): it is not for diagnostic or therapeutic procedures and cannot be purchased by patients.

Isotype Control:

[Bulk Human IgG1 Isotype Control \(IB1\) \[ICH2254\]](#)