



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



## Farletuzumab Biosimilar – Research Grade

### Description

## Farletuzumab Biosimilar – Research Grade

#### Product Benefits:

ichorbio's Farletuzumab Biosimilar for Research Use Only (RUO) and not for therapeutic use.

#### Target:

FOLR1

#### Size:

ichorbio's Farletuzumab Biosimilar is available in the following sizes:  
1mg, 5mg, 10mg, 20mg, 50mg, 100mg

#### Isotype:

Human IgG1(E356D/M358L)-Kappa

#### Host:

CHO Cells

#### Species Reactivity:

Human

#### Background:

Uses the same sequence as the therapeutic antibody Farletuzumab

---

### Specificity:

Detects human FOLR1. This non-therapeutic antibody uses the same sequence as the therapeutic antibody Farletuzumab.

### Purification Method:

This monoclonal antibody was purified using multi-step affinity chromatography methods such as Protein A or G depending on the species and isotype.

### Concentration:

5mg/ml

### Formulation:

Sterile, preservative-free, solution in PBS. BSA and Azide free.

### Purity:

>95% by SDS-PAGE

### Endotoxin:

? 1.0 EU/mg as determined by the LAL method

### Storage:

This antibody is stable for at least 4 weeks when stored at 2-8°C. For long term storage, aliquot in working volumes without diluting and store at – 20°C or -80°C. Avoid repeated freeze thaw cycles.

### Applications:

Functional Assays

### Application Notes:

This biosimilar is strictly for research use only (RUO): it is not for diagnostic or therapeutic procedures and cannot be purchased by patients.