



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Isotype Control IgG2a

Fluorochrome	Reference	Size (µg)
PURE	ICIGG2APU-100	100
FITC	ICIGG2AF-100	100
PE	ICIGG2APE-50	50
APC	ICIGG2AA-50	50
PerCP	ICIGG2APP-100	100
Biotin	ICIGG2AB-100	100
CF- Blue	ICIGG2ACFB-100	100
APC-C750	ICIGG2AAC750-50	50
PerCP-Cyane5.5	ICIGG2APP5.5-100	100

### PRODUCT DESCRIPTION

**Clone:** B12/8

**Isotype:** IgG2a

**Tested application:** flow cytometry

**Host:** Mouse

**Species reactivity:** Human, mouse and rat

**Storage instruction:** store in the dark at 2-8 °C

**Storage buffer:** aqueous buffered solution containing protein stabilizer and 0.09% sodium azide (NaN<sub>3</sub>).

**Recommended usage:** Immunostep's isotype control IgG2a, clone B12/8, is a monoclonal antibody produced against a synthetic hapten, which is normally not present in humans or animals using in flow cytometry. This reagent is effective for direct immunofluorescence staining of human tissue for flow cytometric analysis using  $\leq 1 \mu\text{g}/10^6$  cells.

**Presentation:** liquid

**Source:** Supernatant proceeding from an *in vitro* cell culture of a cell hybridoma.

**Purification:** Affinity chromatography.

### ANTIGEN DETAILS

**Large description:** Today, all techniques used for both diagnosis and research must have negative controls. In Flow Cytometry, these are known as Isotypic Controls. With these it is possible to ensure that our results are reliable.<sup>(1-3)</sup>

The antibody specificity of the mouse IgG2a produced by clone B12/8 was chosen as an isotype control after screening on a variety of live mouse, rat, and human tissues.

Please, refer to [www.immunostep.com](http://www.immunostep.com) technical support for more information.

### WARRANTY

Warranted only to conform to the quantity and contents stated on the label or in the product labelling at the time of delivery to the customer. Immunostep disclaims hereby other warranties.

Immunostep's sole liability is limited to either the replacement of the products or refund of the purchase price.

### REFERENCES

1. Abo Al Hayja M, Eklund A, Grunewald J, Wahlstrom J. Reduced expression of peroxisome proliferator-activated receptor alpha in BAL and blood T cells of non-lofgren's sarcoidosis patients. *J Inflamm (Lond)*;12:28.
2. Jun E, Kim S, Kim JH, Cha K, So IS, Son HN, et al. Design of a multicomponent peptide-woven nanocomplex for delivery of siRNA. *PLoS One*;10(2):e0118310.
3. Tashima Y, Stanley P. Antibodies that detect O-linked beta-D-N-acetylglucosamine on the extracellular domain of cell surface glycoproteins. *J Biol Chem* Apr 18;289(16):11132-42.

### MANUFACTURED BY



#### Immunostep S.L

Avda. Universidad de Coimbra, s/n  
Cancer Research Center (CIC)  
Campus Miguel de Unamuno  
37007 Salamanca (Spain)  
Tel. (+34) 923 294 827  
[www.immunostep.com](http://www.immunostep.com)