



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## NH2-PEG6-Boc

Cat. No.:	HY-130486
CAS No.:	1286281-32-0
Molecular Formula:	C <sub>19</sub> H <sub>39</sub> NO <sub>8</sub>
Molecular Weight:	409.51
Target:	PROTAC Linkers; ADC Linker
Pathway:	PROTAC; Antibody-drug Conjugate/ADC Related
Storage:	4°C, protect from light * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)



### BIOLOGICAL ACTIVITY

<b>Description</b>	NH2-PEG6-Boc is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs <sup>[1]</sup> . NH2-PEG6-Boc is also a non-cleavable 6 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs) <sup>[2]</sup> .		
<b>IC<sub>50</sub> &amp; Target</b>	PEGs	Alkyl/ether	Non-cleavable Linker
<b>In Vitro</b>	PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins <sup>[1]</sup> . ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker <sup>[2]</sup> . MCE has not independently confirmed the accuracy of these methods. They are for reference only.		

### REFERENCES

- [1]. Foley CA, et al. Assessing the Cell Permeability of Bivalent Chemical Degraders Using the Chloroalkane Penetration Assay. ACS Chem Biol. 2020 Jan 17;15(1):290-295.
- [2]. Kenneth John DIRICO, et al. Spliceostatin analogs. WO2014068443A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA