



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Amino-PEG4-CH₂COOH

Cat. No.:	HY-130524
CAS No.:	195071-49-9
Molecular Formula:	C ₁₀ H ₂₁ NO ₆
Molecular Weight:	251.28
Target:	PROTAC Linkers; ADC Linker
Pathway:	PROTAC; Antibody-drug Conjugate/ADC Related
Storage:	4°C, protect from light * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)



SOLVENT & SOLUBILITY

In Vitro

DMSO : 100 mg/mL (397.96 mM; Need ultrasonic)

Concentration	Solvent	Mass		
		1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM	3.9796 mL	19.8981 mL	39.7962 mL
	5 mM	0.7959 mL	3.9796 mL	7.9592 mL
	10 mM	0.3980 mL	1.9898 mL	3.9796 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description

Amino-PEG4-CH₂COOH is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs^[1]. Amino-PEG4-CH₂COOH is also a non-cleavable 4 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)^[2].

IC₅₀ & Target

PEGs Non-cleavable Linker

In Vitro

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins^[1]. ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker^[2]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Nan JI, et al. Merck degraders and uses thereof. WO2020010210A1.

[2]. Maggie Lu, et al. Linker-drug and antibody-drug conjugate (adc) employing the same. EP3335734A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA