



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mal-Gly-PAB-Exatecan-D-glucuronic acid

Cat. No.:	HY-153179
CAS No.:	2763252-25-9
Molecular Formula:	C ₄₇ H ₄₅ FN ₆ O ₁₇
Molecular Weight:	984.89
Target:	Drug-Linker Conjugates for ADC
Pathway:	Antibody-drug Conjugate/ADC Related
Storage:	4°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)



SOLVENT & SOLUBILITY

In Vitro

DMSO : 100 mg/mL (101.53 mM; Need ultrasonic)

Concentration	Mass		
	1 mg	5 mg	10 mg
1 mM	1.0153 mL	5.0767 mL	10.1534 mL
5 mM	0.2031 mL	1.0153 mL	2.0307 mL
10 mM	0.1015 mL	0.5077 mL	1.0153 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description

Mal-Gly-PAB-Exatecan-D-glucuronic acid can be used for Antibody-Drug Conjugates (ADCs) synthesis. Mal-Gly-PAB-Exatecan-D-glucuronic acid can be used for research of tumor research and diagnosis^[1].

REFERENCES

[1]. Anderl Jan, et al. Anti-CEACAM5 antibodies and conjugates and uses thereof. World Intellectual Property Organization, WO2022048883 A1. 2022-03-10.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA