



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Cross-linked dextran G 200

---

Cat. No.: HY-156045C  
CAS No.: 9041-36-5  
Target: Others  
Pathway: Others  
Storage: 4°C, sealed storage, away from moisture

Cross-linked dextran G 200

### BIOLOGICAL ACTIVITY

#### Description

Cross-linked dextran G 200 is a hydrophilic gel that can be used as a gel filter filler (Particle size range: 40-120 µm; Sphere protein separation range: 5k-600k Da)<sup>[1]</sup>.

### REFERENCES

[1]. Ogston AG, et al. The thermodynamics of interaction between Sephadex and penetrating solutes. *Biochem J.* 1970 Jan;116(2):171-5.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: [tech@MedChemExpress.com](mailto:tech@MedChemExpress.com)

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA