



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Flutax-2

Cat. No.:	HY-162161
CAS No.:	301844-13-3
Molecular Formula:	C ₇₁ H ₆₄ F ₂ N ₂ O ₂₁
Molecular Weight:	1319.27
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



BIOLOGICAL ACTIVITY

Description

Flutax-2 is an active fluorescent derivative of Paclitaxel, binds to $\alpha\beta$ -tubulin dimer polymerized. Flutax-2 can be used for imaging microtubules in live cells, isolated cytoskeletons and parasite (Ex/Em=496/526 nm)^[1].

REFERENCES

[1]. Sheila B Lecke, et al. *Trichomonas vaginalis*: microtubule cytoskeleton distribution using fluorescent taxoid. *Exp Parasitol*. 2002 Oct;102(2):113-6.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA