



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BODIPY-FL

Cat. No.:	HY-43520
CAS No.:	165599-63-3
Molecular Formula:	C ₁₄ H ₁₅ BF ₂ N ₂ O ₂
Molecular Weight:	292.09
Target:	Fluorescent Dye
Pathway:	Others
Storage:	4°C, sealed storage, away from moisture and light * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)



SOLVENT & SOLUBILITY

In Vitro

DMSO : 100 mg/mL (342.36 mM; Need ultrasonic)

Concentration	Solvent	Mass		
		1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM	3.4236 mL	17.1180 mL	34.2360 mL
	5 mM	0.6847 mL	3.4236 mL	6.8472 mL
	10 mM	0.3424 mL	1.7118 mL	3.4236 mL

Please refer to the solubility information to select the appropriate solvent.

BIOLOGICAL ACTIVITY

Description

BODIPY-FL is a potent fluorescent dye. BODIPY-FL can be used to label probe or primer for fluorescent quenching-based quantitative detection of specific DNA/RNA. BODIPY-FL-labeled monoterpenoid can be used to examine the features of a broad spectrum of Gram-positive and Gram-negative bacteria and pathogenic fungi as well^{[1][2][3]}.

REFERENCES

- [1]. David Dulin, et al. Reduced photobleaching of BODIPY-FL. *Physics Procedia*. 2010, February 28, 1563-1567.
- [2]. Kurata S, et al. Fluorescent quenching-based quantitative detection of specific DNA/RNA using a BODIPY((R)) FL-labeled probe or primer. *Nucleic Acids Res*. 2001 Mar 15;29(6):E34.
- [3]. Guseva GB, et al. Design, Spectral Characteristics, and Possibilities for Practical Application of BODIPY FL-Labeled Monoterpenoid. *ACS Appl Bio Mater*. 2021 Aug 16;4(8):6227-6235.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA