



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## 6-Iodoacetamidofluorescein

Cat. No.:	HY-D0122
CAS No.:	73264-12-7
Molecular Formula:	C <sub>22</sub> H <sub>14</sub> INO <sub>6</sub>
Molecular Weight:	515.25
Target:	Fluorescent Dye
Pathway:	Others
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.



### BIOLOGICAL ACTIVITY

#### Description

6-Iodoacetamidofluorescein, a sulphhydryl-specific fluorescent dye, and can be used to selectively label the -SH groups of nuclear matrix polypeptides and proteins<sup>[1][2]</sup>.

### REFERENCES

[1]. Martelli AM, et, al. Probing the sulfhydryl groups of nuclear matrix proteins with 6-iodoacetamidofluorescein. *Cell Biol Int Rep.* 1990 May;14(5):409-18.

[2]. Cocco L, et, al. Changes in nucleosome structure and histone H3 accessibility. Iodoacetamidofluorescein labelling after treatment with phosphatidylserine vesicles. *Exp Cell Res.* 1986 Oct;166(2):465-74.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA