



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## 18-Crown-6-ether

Cat. No.:	HY-D0180
CAS No.:	17455-13-9
Molecular Formula:	C <sub>12</sub> H <sub>24</sub> O <sub>6</sub>
Molecular Weight:	264.32
Target:	Biochemical Assay Reagents
Pathway:	Others
Storage:	4°C, sealed storage, away from moisture * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)



### BIOLOGICAL ACTIVITY

#### Description

18-Crown-6-ether (18C6; 1,4,7,10,13,16-Hexaoxacyclooctadecane) is an organic compound that can be used to prepare stable microcapsule responsive layers for further assembly into bilayer microcapsules. For example, 18-Crown-6-ether is used to prepare the response layer and is coated with a G-quadruplex cross-linked hydrogel layer stabilized by K<sup>+</sup>; when Mg<sup>2+</sup> ions are present, 18-Crown-6-ether and K<sup>+</sup> ions can respectively dissociate and lock with the G-quadruplex cross-linked layer, thereby achieving switchable controlled release of the load<sup>[1]</sup>.

### REFERENCES

[1]. Lilienthal S, et al. Single and Bilayer Polyacrylamide Hydrogel-Based Microcapsules for the Triggered Release of Loads, Logic Gate Operations, and Intercommunication between Microcapsules. ACS Appl Mater Interfaces. 2020 Jul 15;12(28):31124-31136.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA