



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Bovine Serum Albumin (Cell culture, Endotoxin<0.01EU/mg)

Cat. No.:	HY-D0842A		
CAS No.:	9048-46-8		
Target:	Biochemical Assay Reagents		
Pathway:	Others		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month

Bovine Serum Albumin (Cell culture, Endotoxin<0.01EU/mg)

### SOLVENT & SOLUBILITY

**In Vitro** H<sub>2</sub>O : 200 mg/mL (Need ultrasonic)

### BIOLOGICAL ACTIVITY

**Description** Bovine Serum Albumin (Cell culture, Endotoxin<0.01EU/mg) (BSA) is a 583-residue protein consisting of three homologous all- $\alpha$  domains, organized in a heart-shaped structure. BSA is a globular protein that is used in numerous biochemical applications.

### REFERENCES

[1]. Holm NK, et al. Aggregation and fibrillation of bovine serum albumin. *Biochim Biophys Acta*. 2007 Sep;1774(9):1128-38.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA