



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

FITC-GW3965

Cat. No.:	HY-D1734
CAS No.:	2374144-23-5
Molecular Formula:	C ₅₉ H ₅₆ N ₄ O ₉ S
Molecular Weight:	997.16
Target:	LXR
Pathway:	Metabolic Enzyme/Protease; Vitamin D Related/Nuclear Receptor
Storage:	4°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)



BIOLOGICAL ACTIVITY

Description

FITC-GW3965 is a fluorescence-labelled liver X receptor β (LXR β) agonist GW3965 (HY-10627). FITC-GW3965 is a tracer, that can be designed by replacing the trifluoromethyl of GW3965 with an amide to link the FITC. FITC-GW3965 can be used to study the function of LXR β ^[1].

IC₅₀ & Target

LXR β ^[1]

REFERENCES

[1]. Zhang Z, et al. Identify liver X receptor β modulator building blocks by developing a fluorescence polarization-based competition assay. Eur J Med Chem. 2019 Sep 15;178:458-467.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA