



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Formimino-L-glutamic acid Transferase, porcine liver

Cat. No.:	HY-E70399	
CAS No.:	9032-83-1	
Target:	Biochemical Assay Reagents	
Pathway:	Others	Formimino-L-glutamic acid Transferase, porcine liver
Storage:	Please store the product under the recommended conditions in the Certificate of Analysis.	

### BIOLOGICAL ACTIVITY

Description	Formimino-L-glutamic acid Transferase can convert aminomethyleneglutamic acid FIGLU and tetrahydrofolate THF into L-glutamic acid and 5-formimino-THF <sup>[1]</sup> .
-------------	--

### REFERENCES

[1]. Lv G, et al. Marker-assisted pyramiding of  $\gamma$ -tocopherol methyltransferase and glutamate formiminotransferase genes for development of biofortified sweet corn hybrids. PeerJ. 2022 Jul 6;10:e13629.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA