



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Phytohemagglutinin P

Cat. No.:	HY-N7038A
Target:	Others
Pathway:	Others
Storage:	4°C, protect from light * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)

Phytohemagglutinin P

SOLVENT & SOLUBILITY

In Vitro	H ₂ O : 10 mg/mL (Need ultrasonic and warming)
In Vivo	1. Add each solvent one by one: PBS Solubility: 25 mg/mL (Infinity mM); Clear solution; Need ultrasonic

BIOLOGICAL ACTIVITY

Description	Phytohemagglutinin P (PHA-P) is a mitogen known to selectively stimulate cells of hematogenous or lymphoid monocytic origin ^[1] .
In Vitro	Phytohemagglutinin-P (10 µg/ml; 24-48hous)-stimulated decidual CD8 ⁺ T lymphocytes produces high levels of both interferon gamma and interleukin (IL) 8, and low levels of granulocyte-macrophage colony-stimulating factor (CSF2), IL1B, IL2, IL6, IL10, IL12, and tumor necrosis factor; these did not vary with gestational age ^[2] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.
In Vivo	Phytohemagglutinin-P (25 and 50 mg/kg; i.p.; daily for 15 days; male Sprague-Dawley rats) shows a decrease in longitudinal bone growth, in the number of chondrocytes, in the thickness of cartilage plate, in the metaphyseal mass of hard tissue, in the percentage of calcified cartilage core, and in the number of osteoblasts per mm of bone surface. Additionally, PHA-P increases the number of osteoclasts, the number of labeled osteoclastic nuclei, and the average number of nuclei per osteoclast ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Scaife PJ, et al. Effector activity of decidual CD8⁺ T lymphocytes in early human pregnancy. Biol Reprod. 2006;75(4):562-567.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA