



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## VTSEGAGLQK-<sup>13</sup>C<sub>6</sub>,<sup>15</sup>N<sub>2</sub>

<b>Cat. No.:</b>	HY-P5201S
<b>Molecular Formula:</b>	C <sub>46</sub> <sup>13</sup> C <sub>6</sub> H <sub>91</sub> N <sub>13</sub> <sup>15</sup> N <sub>2</sub> O <sub>19</sub>
<b>Molecular Weight:</b>	1238.31
<b>Target:</b>	Isotope-Labeled Compounds
<b>Pathway:</b>	Others
<b>Storage:</b>	Please store the product under the recommended conditions in the Certificate of Analysis.



### BIOLOGICAL ACTIVITY

<b>Description</b>	VTSEGAGLQK- <sup>13</sup> C <sub>6</sub> , <sup>15</sup> N <sub>2</sub> is the <sup>13</sup> C and <sup>15</sup> N labeled VTSEGAGLQK. VTSEGAGLQK is an amino acid sequence in the C-terminal region of recombinant human alpha-acid glucosidase (rhGAA). VTSEGAGLQK can be used to bind anti-drug antibodies (ADA) in plasma and quantitatively analyze the therapeutic effect <sup>[1][2]</sup> .
<b>In Vitro</b>	Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs <sup>[1]</sup> . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### REFERENCES

- [1]. Russak EM, et al. Impact of Deuterium Substitution on the Pharmacokinetics of Pharmaceuticals. *Ann Pharmacother.* 2019 Feb;53(2):211-217.
- [2]. Bronsema KJ, et al. Absolute quantification of the total and antidrug antibody-bound concentrations of recombinant human  $\alpha$ -glucosidase in human plasma using protein G extraction and LC-MS/MS. *Anal Chem.* 2015 Apr 21;87(8):4394-401.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA