



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

GLYATL1 Human Pre-designed siRNA Set A

| | | | |
|-----------|-------------------------------|-------|----------|
| Cat. No.: | HY-RS16225 | | |
| Target: | Small Interfering RNA (siRNA) | | |
| Pathway: | Epigenetics | | |
| Storage: | Powder | -20°C | 3 years |
| | In solvent | -80°C | 6 months |
| | | -20°C | 1 month |

GLYATL1 Human Pre-designed siRNA Set A

BIOLOGICAL ACTIVITY

Description

GLYATL1 Human Pre-designed siRNA Set A contains three designed siRNAs for GLYATL1 gene (Human), as well as a negative control, a positive control, and a FAM-labeled negative control.

Components

GLYATL1 siRNA-1: 5 nmol (HPLC)

GLYATL1 siRNA-2: 5 nmol (HPLC)

GLYATL1 siRNA-3: 5 nmol (HPLC)

siRNA Negative Control: 5 nmol (HPLC)

FAM-labeled siRNA Negative Control: 5 nmol (HPLC)

GAPDH siRNA Positive Control: 5 nmol (HPLC)

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA