



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## 4-Bromomethyl-6,7-dimethoxycoumarin

Cat. No.:	HY-W269179
CAS No.:	88404-25-5
Molecular Formula:	C <sub>12</sub> H <sub>11</sub> BrO <sub>4</sub>
Molecular Weight:	299.12
Target:	Fluorescent Dye
Pathway:	Others
Storage:	4°C, protect from light * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)



### SOLVENT & SOLUBILITY

#### In Vitro

DMSO : 10 mg/mL (33.43 mM; ultrasonic and warming and heat to 70°C)

Concentration	Solvent	Mass		
		1 mg	5 mg	10 mg
Preparing Stock Solutions	1 mM	3.3431 mL	16.7157 mL	33.4314 mL
	5 mM	0.6686 mL	3.3431 mL	6.6863 mL
	10 mM	0.3343 mL	1.6716 mL	3.3431 mL

Please refer to the solubility information to select the appropriate solvent.

### BIOLOGICAL ACTIVITY

#### Description

4-Bromomethyl-6,7-dimethoxycoumarin is a fluorescent label for carboxylic acids in chromatographic detection<sup>[1]</sup>.

### REFERENCES

[1]. Amet Y, et al. Simultaneous radiometric and fluorimetric detection of lauric acid metabolites using high-performance liquid chromatography following esterification with 4-bromomethyl-6,7-dimethoxycoumarin in human and rat liver microsomes. J Chromatogr B Biomed Appl. 1996 Jun 7;681(2):233-9.

[2]. R.Farinotti, et al. 4-Bromomethyl-6,7-dimethoxycoumarin as a fluorescent label for carboxylic acids in chromatographic detection. Journal of Chromatography A. Volume 269, 1983, Pages 81-90.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA