



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Isethionic acid

Cat. No.:	HY-Y0095
CAS No.:	107-36-8
Molecular Formula:	C ₂ H ₆ O ₄ S
Molecular Weight:	126.13
Target:	Endogenous Metabolite
Pathway:	Metabolic Enzyme/Protease
Storage:	Solution, -20°C, 2 years



SOLVENT & SOLUBILITY

In Vitro	DMSO : 100 mg/mL (792.83 mM; Need ultrasonic)
In Vivo	<ol style="list-style-type: none"> Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (19.82 mM); Clear solution Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (19.82 mM); Clear solution Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (19.82 mM); Clear solution

BIOLOGICAL ACTIVITY

Description	Isethionic acid (2-Hydroxyethanesulfonic acid) is an organosulfur compound. Isethionic acid is widely distributed in animal species and in a few red algal species. Isethionic acid can be used as an anionic detergent and has anti-settlement activity against <i>Balanus amphitrite</i> ^{[1][2]} .
IC ₅₀ & Target	Human Endogenous Metabolite

REFERENCES

- [1]. Hellio C, et al. Isethionic acid and floridoside isolated from the red alga, *Grateloupia turuturu*, inhibit settlement of *Balanus amphitrite* cyprid larvae. *Biofouling*. 2004 Jun;20(3):139-45.
- [2]. Tupker RA, et al. Irritancy ranking of anionic detergents using one-time occlusive, repeated occlusive and repeated open tests. *Contact Dermatitis*. 1999 Jun;40(6):316-22.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA