



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7


mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Farnesol-d₆

Cat. No.:	HY-Y0248AS		
CAS No.:	166447-71-8		
Molecular Formula:	C ₁₅ H ₂₀ D ₆ O		
Molecular Weight:	228.4		
Target:	Bacterial; Endogenous Metabolite; Antibiotic; Isotope-Labeled Compounds		
Pathway:	Anti-infection; Metabolic Enzyme/Protease; Others		
Storage:	Pure form	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month



BIOLOGICAL ACTIVITY

Description	Farnesol-d ₆ is deuterium labeled Farnesol. Farnesol is a sesquiterpene alcohol that modulates cell-to-cell communication in <i>Candida albicans</i> , and has the activity in inhibiting bacteria.
In Vitro	Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

- [1]. Russak EM, et al. Impact of Deuterium Substitution on the Pharmacokinetics of Pharmaceuticals. *Ann Pharmacother.* 2019;53(2):211-217.
- [2]. Cordeiro Rde A, et al. Farnesol inhibits in vitro growth of the *Cryptococcus neoformans* species complex with no significant changes in virulence-related exoenzymes. *Vet Microbiol.* 2012 Oct 12;159(3-4):375-80.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA