

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Proteins

Screening Libraries

Product Data Sheet



Cat. No.: HY-Y1069 CAS No.: 97-67-6 Molecular Formula: $C_4H_6O_5$ Molecular Weight: 134

Target: **Endogenous Metabolite** Pathway: Metabolic Enzyme/Protease

Storage: Powder -20°C 4°C

2 years In solvent -80°C 6 months

> -20°C 1 month

3 years

SOLVENT & SOLUBILITY

In Vitro

H₂O: 100 mg/mL (746.27 mM; Need ultrasonic) DMSO: 100 mg/mL (746.27 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	7.4627 mL	37.3134 mL	74.6269 mL
	5 mM	1.4925 mL	7.4627 mL	14.9254 mL
	10 mM	0.7463 mL	3.7313 mL	7.4627 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (18.66 mM); Clear solution
- 2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (18.66 mM); Clear solution
- 3. Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (18.66 mM); Clear solution

BIOLOGICAL ACTIVITY

Description	(S)-Malic acid ((S)-2-Hydroxysuccinic acid) is a dicarboxylic acid in naturally occurring form, contributes to the pleasantly sour taste of fruits and is used as a food additive.	
IC ₅₀ & Target	Human Endogenous Metabolite	Human Endogenous Metabolite

In Vitro

It is showed that ME is essential for (S)-2-Hydroxysuccinic acid (L-malic acid) utilization in L. casei. Furthermore, deletion of either the gene encoding the histidine kinase or the response regulator of the TC system resulted in the loss of the ability to grow on (S)-2-Hydroxysuccinic acid, thus indicating that the cognate TC system regulates and is essential for the expression of ME. Transcriptional analyses shows that expression of maeE is induced in the presence of (S)-2-Hydroxysuccinic acid and repressed by glucose, whereas TC system expression is induced by (S)-2-Hydroxysuccinic acid and is not repressed by glucose^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

REFERENCES

[1]. Landete JM, et al. Requirement of the Lactobacillus casei MaeKR two-component system for L-malic acid utilization via a malic enzyme pathway. Appl Environ Microbiol. 2010 Jan;76(1):84-95.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898 Fax: 609-228-5909 E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA