



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



## Pooled mouse ascites without immunoglobulins

Catalogue number: **NOR/AA**

Product Type	Serum
Units	0.5 ml
Species reactivity	Mouse

### Background

Ascites obtained from mice bearing different hybridomas, pooled and deleted of all immunoglobulins. It is used as a negative control for ELISA and for indirect immunofluorescence and immunoperoxidase methods, replacing the specific ascites in the test arrangement. The control ascites is used in a dilution comparable to that of the specific ascites.

### Product

Delipidated, heat inactivated, lyophilized, stable pooled ascites. No preservative added.

### Applications

The lyophilized ascites is shipped at ambient temperature. After arrival store at +4°C. Prolonged storage may be at -20°C. After opening the vial, the lyophilized content is reconstituted by adding 1 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted ascites may be stored frozen.

### Specificity

Pooled normal Mouse Serum

### Storage

Lyophilized at +4° C--at least 10 years. Reconstituted at or below -20° C--3-5 years. Reconstituted at +4° C--7 days.

### Caution

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

## Distributors

For Purchasing Information, please contact your local distributor

[Find Distributor](#)

