



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Human IgA1 kappa, polymer

Catalogue number: **P 430**

Product Type	Proteins & Peptides
Units	0.5 mg
Species reactivity	Human

Background

Confirmed by immunoelectrophoresis and double radial immunodiffusion (Ouchterlony). Not less than 98% as determined by agar-gel electrophoresis, quantitative densitometry immunoelectrophoresis and double radial immunodiffusion (Ouchterlony). This IgA1 consists of polymers and will contain J chain. The product contains a trace of non-Ig protein. Purified human IgA1 kappa may be used as a reference antigen, calibrator, coating protein and blocking agent in a variety of immunoassays including immunodiffusion, immunoelectrophoresis, ELISA, Western blotting, dot-immunobinding assays (DIBA), haemagglutination, and cell-binding assays.

Product

Lyophilized, stable white powder No preservative added

Applications

Specificity

Chromatographically purified Human IgA1 kappa, polymer.

Storage

This product is shipped lyophilized at ambient temperature. After arrival store at +4°C. Prolonged storage may be at -20°C. The lyophilized content is reconstituted by adding 0.5 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted solution may be stored frozen. Working dilutions may be prepared by adding the required amount of phosphate buffered saline (PBS), should not be refrozen and preferably used within 24 hours. Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation and will not affect the performance of the product. Lyophilized at +4°C--at least 10 years. Reconstituted at or below -20°C--3-5 years. Reconstituted at +4°C--7 days.

Caution

Distributors

For Purchasing Information, please contact your local distributor

[Find Distributor](#)

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUbio accepts no liability for any inaccuracies or omissions in this information.

[Home](#) | [Company Profile](#) | [Catalogue](#) | [Distributors](#) | [Contact](#) |

Content: Nordic-MUbio BV - Copyright © 2015