

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



## Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

## SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





www.rockland.com tech@rockland.com +1 484.791.3823

### Datasheet for 001-G70-MF0

# **Ultrapure Collagen I for plate coating - Bovine Placenta**

### **Overview**

Description:	Ultrapure Collagen I for plate coating - Bovine Placenta - 001-G70-MF0
Item No.:	001-G70-MF0
Size:	8 mL
Origin:	Bovine

## **Product Details**

Background:	Bovine Collagen Type I is one of a family of proteins found particularly in the flesh and connective tissues of mammals (approximately one-third of the body's total protein). Over two dozen types of collagen have been described; Type I is the most abundant form in the body. Collagen Type I is found in scar tissue, tendons, the skin, arterial walls, the corneas, muscles, cartilage, and in the organic parts of bones and teeth. Bovine Collagen Type I is ideal for investigators involved in extracellular matrix proteins and osteoporosis.
Synonyms:	Type I collagen, collagen 1, bovine collagen, Alpha 1 type I collagen, Alpha 2 type I collagen, COL1A1, COL1A2, Collagen alpha-1(I) chain, Collagen alpha-2 (I) chain, plate coating solution
Species of Origin:	Bovine

## **Target Details**

Gene Name:	COL1A1
Purity/Specificity:	Ultrapure Bovine Collagen Type I Solution is purified from bovine placenta and shown to be 99% pure by gel electrophoresis and immunoassays.
Relevant Links:	Collagen_Top_Culture-Protocol

# **Application Details**

Application Note:	Ultrapure Bovine Collagen Type I solution are suitable for plate coating. Enough to coat 10 (ten) T75 plates.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

www.rockland.com Page 1 of 3





www.rockland.com tech@rockland.com +1 484.791.3823

#### **Formulation**

Physical State:	Liquid (sterile filtered)
Concentration:	1.0 mg/mL by dry weight
Buffer:	0.013M HCl
Preservative:	None
Stabilizer:	None

## **Shipping & Handling**

<b>Shipping Condition:</b>	Wet Ice
Storage Condition:	Store vial at 4° C prior to opening. This product is stable at 4° C as an undiluted liquid. Dilute only prior to immediate use. Keep under sterile conditions.
Expiration:	Expiration date is six (6) months from date of receipt.

### **Disclaimer**

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

www.rockland.com Page 2 of 3





www.rockland.com tech@rockland.com +1 484.791.3823

www.rockland.com Page 3 of 3