



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet for TMS-15-50**15/50 Tube TrueBlot® Magnetic Separator****Overview**

Description:	15/50 Tube TrueBlot® Magnetic Separator - TMS-15-50
Item No.:	TMS-15-50
Size:	1 Each

Product Details

Background:	A magnetic separation unit designed to hold, and separate magnetic particles to the side of the tubes, up to five 15 mL and three 50 mL tubes. The unit contains five bar magnets embedded in a plastic housing.
Synonyms:	holder, rack, magnets, magnet

Application Details

Application Note:	Magnetic separator unit can be used to hold, and separate magnetic particles.
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.


Formulation

Concentration:	1 Each
-----------------------	--------

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store @ room temperature (18° - 26° C)

Images

**Separator**

This magnetic separator is designed to hold, and separate magnetic particles to the side of the tubes, up to five 15 mL and three 50 mL tubes.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.



Order online now!

www.rockland.com
tech@rockland.com
+1 484.791.3823