

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





www.rockland.com tech@rockland.com +1 484.791.3823

Datasheet for WM3702-02-1000

WM3702 Purified Genomic DNA

Overview

Description:	WM3702 Purified Genomic DNA - WM3702-02-1000
Item No.:	WM3702-02-1000
Size:	10 μg
Origin:	Human

Product Details

Background: DNA was prepared from cell line WM3702. WM3702 is a metastatic human melanoma cell line

that displays pigmented, thin elongated morphology in culture. This cell line was established from a right upper flank of back in a 58-year-old male with necrotic fungating melanoma. The subject displayed metastatic melanoma forming multiple subcutaneous nodules with overlying skin ulcerations. This cell line contains a Q61L mutation at position 61 in the N-RAS gene. The Q61L mutation results in an amino acid substitution at position 61 in NRAS, from a glutamine (Q)

to a leucine (L). WM3702 cells match with WM3758 and WM3682. These cells produce $\,$

xenograft tumors when injected into immunocompromised mice.

Synonyms: Melanoma, patient derived tumor, tumor models, skin cancer, xenograft

Species of Origin: Human

Target Details

Purity/Specificity: Genomic DNA was purified from cells using genomic DNA preparaton kit according to

manufacturers instruction. DNA was diluted to 200 $ng/\mu L$ in TE buffer (0.01 M Tris Chloride, 0.001 M EDTA, pH 7.6). Concentration was determined at A260 using nanodrop ND-1000.

0.001 M LDTA, ph 7.0). Concentration was determined at A200 using nanourop ND-1000

Relevant Links: • Cell Line EULA

Application Details

Application Note: Purified Genomic DNA is suitable for a number of molecular biology applications including but

not limited to preparation of genomic libraries, PCR templates, DNA sequencing, DNA

fingerprinting, and mutation analysis.

www.rockland.com Page 1 of 3





www.rockland.com tech@rockland.com +1 484.791.3823

Assay Dilutions: All assays should be optimized by the user. Recommended dilutions (if any) may be

listed below.

Cell Line Data

Cell Line:	Human Melanoma
Product Type:	DNA
Cell Viability:	No
Stage:	Metastasis
BRAF:	WT
CDK4:	WT
C-Kit:	WT
N-RAS:	Q61L
PTEN:	ND
Paired:	Yes

Formulation

Physical State:	Liquid
Concentration:	200 ng / μl by UV absorbance at 260 nm
Buffer:	0.01 M Tris Chloride, 0.001 M EDTA, pH 7.6
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Dry Ice
Storage Condition:	Store vial at -20° C or COLDER. For extended storage, aliquot contents to minimize freeze/thaw cycles.
Expiration:	Expiration date is six (6) months from date of receipt.

Images

www.rockland.com Page 2 of 3





www.rockland.com tech@rockland.com +1 484.791.3823



VialHuman melanoma tumor cell genomic DNA isolated with genomic DNA miniprep kit

Disclaimer

No test method can provide total assurance that the hepatitis B virus, hepatitis C virus, human immunodeficiency virus, or any other infectious agents are absent. Thus, all blood products, including purified proteins derived from human blood sources, should be handled at Biosafety Level 2 as recommended by the CDC\NIH manual entitled Biosafety in Microbiological and Biomedical Laboratories for potentially infectious human serum, blood specimens or proteins derived from same. Source material for the human blood product supplied to your facility has been tested for the detection of HIV antibody, Hepatitis B surface antigen, antibody to Hepatitis C, HIV 1 antigen(s), antibody to HTLV - I/II, and syphilis by FDA guidelines. All units were found to be non-reactive/negative for these tests. All human blood source material is collected in FDA licensed centers and is tested with FDA approved test kits.

Cell Line Limited Use License Required. THIS PRODUCT IS SUBJECT TO AN END-USER LICENSE AGREEMENT (EULA). BY ACCEPTING THIS PRODUCT, RECIPIENT AGREES TO BE BOUND BY THE TERMS OF USE SET FORTH BELOW and SET FORTH IN THE EULA. THIS PRODUCT IS FOR IN VITRO RESEARCH USE ONLY. THERAPEUTIC, DIAGNOSTIC, OR VETERINARY USE IS PROHIBITED. This product may not be resold or transferred by the recipient and may be used only by the recipient, in the recipient's facility and only for research use and other uses specifically permitted by the EULA. No other commercial use is allowed. "Commercial Use" means any and all uses of this product by recipient or others for monetary or other consideration, including providing services, supplying information or data to unaffiliated third parties, and resale or transfer of this product for any use. Recipient has no right to modify, derivatize, genetically engineer or otherwise create variations of this product or associated cells or cell lines. ROCKLAND AND WISTAR MAKE NO REPRESENTATIONS AND EXTEND NO WARRANTIES OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING AS TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE PRODUCTS WILL NOT INFRINGE ANY PATENT, COPYRIGHT, TRADEMARK, OR OTHER PROPRIETARY RIGHTS. The terms set forth herein and in the EULA shall be governed by the laws of the Commonwealth of Pennsylvania, USA. To obtain a COMMERCIAL USE license for this product, please contact Rockland Immunochemicals, Inc. Please contact a technical service representative for more information. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

www.rockland.com Page 3 of 3