



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# GlyR $\beta$ (m): 293T Lysate: sc-120533

## BACKGROUND

In the central nervous system (CNS), glycine-mediated inhibitory neurotransmission is essential to voluntary motor control and reflex responses. Glycine binds to glycine receptors (GlyR) in the postsynaptic neuronal membranes. GlyR,  $\gamma$ -aminobutyric acid, serotonin and acetylcholine comprise an evolutionally conserved superfamily of ligand-gated ion channels. The pentameric subunit structure of GlyR consists of two types of glycosylated membrane proteins,  $\alpha$ 1 through  $\alpha$ 4 and  $\beta$ , and an associated peripheral membrane protein, which combine to form a chloride-selective ion channel. In humans, the composition of the pentamer changes from  $\alpha$ 2 subunits in the fetal CNS to  $\alpha$ 1 and  $\beta$  subunits in the adult CNS. Fast potentiation of GlyR by intracellular  $Ca^{2+}$  in the brainstem and midbrain indicate an important role for  $Ca^{2+}$  in modulation of glycinergic synapses.

## REFERENCES

1. Pfeiffer, F., et al. 1981. Solubilisation of the glycine receptor from rat spinal cord. *Brain Res.* 226: 273-279.
2. Pfeiffer, F., et al. 1982. Purification by affinity chromatography of the glycine receptor of rat spinal cord. *J. Biol. Chem.* 257: 9389-9393.
3. Genningsloh, G., et al. 1987. The strychnine-binding subunit of the glycine receptor shows homology with nicotinic acetylcholine receptors. *Nature* 328: 215-220.
4. Schofield, P.R., et al. 1987. Sequence and functional expression of the GABA<sub>A</sub> receptor shows a ligand-gated receptor super-family. *Nature* 328: 221-227.
5. Langosch, D., et al. 1988. Conserved quaternary structure of ligand-gated ion channels: the postsynaptic glycine receptor is a pentamer. *Proc. Natl. Acad. Sci. USA* 85: 7394-7398.
6. Hoch, W., et al. 1989. Primary cultures of mouse spinal cord express the neonatal isoform of the inhibitory glycine receptor. *Neuron* 3: 339-348.
7. Genningsloh, G., et al. 1990. Alpha subunit variants of the human glycine receptor: primary structures, functional expression and chromosomal location of corresponding genes. *EMBO J.* 9: 771-776.
8. Wasmuth, J.J. 1992. A comparison of three methods to produce a high resolution physical map of 11 genes on the distal region of the long arm of human chromosome 5: radiation hybrid mapping, pulsed-field gel electrophoresis and fluorescent *in situ* hybridization. *Am. J. Hum. Genet.* 51: A248.
9. Nikolic, Z., et al. 1998. The human glycine receptor subunit  $\alpha$ -3: GLRA3 gene structure, chromosomal localization, and functional characterization of alternative transcripts. *J. Biol. Chem.* 273: 19708-19714.

## CHROMOSOMAL LOCATION

Genetic locus: Grlb (mouse) mapping to 3 E3.

## PRODUCT

GlyR  $\beta$  (m): 293T Lysate represents a lysate of mouse GlyR  $\beta$  transfected 293T cells and is provided as 100  $\mu$ g protein in 200  $\mu$ l SDS-PAGE buffer.

## APPLICATIONS

GlyR  $\beta$  (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive GlyR  $\beta$  antibodies. Recommended use: 10-20  $\mu$ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

## STORAGE

Store at  $-20^{\circ}$  C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.