



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

GRP 94 (m): 293T Lysate: sc-120647

BACKGROUND

Heat shock protein (HSP) molecular chaperones are environmental stress-inducible gene products. The human HSP 90 family includes 17 genes that fall into four classes: HSP90AA, HSP90AB, HSP90B and TRAP. HSP 90 family members guide the normal folding, intracellular disposition and proteolytic turnover of many key regulators of cell growth, differentiation and survival. HSP 90 α , also designated HSP90A, HSP 86 and LPS-associated protein 2 (LAP2), is a cytosolic enhancer of inducible nitric-oxide synthase (iNOS), with chaperone activity that is important for the transcriptional activity of p53. HSP 90 β , also designated HSP90B, HSP 84 and HSPC2, is a cytosolic protein that participates in signaling pathways with PKC ϵ to protect cells from external damage, particularly in heat shock-mediated events. GRP 94, also known as tumor rejection antigen 1 (TRA1), ECGP and GP96, localizes to the ER, is highly expressed in BGC-823 human gastric carcinoma cells and is upregulated in human endothelial cells in response to hypoxia by HIF-1. TRAP1 (TNF receptor-associated protein 1), also designated HSP 75, is a mitochondrial matrix component that plays a role in the induction of apoptosis in response to reactive oxygen species.

REFERENCES

1. Wu, J.M., Xiao, L., Cheng, X.K., Cui, L.X., Wu, N.H. and Shen, Y.F. 2003. PKC ϵ is a unique regulator for HSP 90 β gene in heat shock response. *J. Biol. Chem.* 278: 51143-51149.
2. Whitesell, L. and Lindquist, S.L. 2005. HSP 90 and the chaperoning of cancer. *Nat. Rev. Cancer* 5: 761-772.
3. Cowen, L.E. and Lindquist, S. 2005. HSP 90 potentiates the rapid evolution of new traits: drug resistance in diverse fungi. *Science* 309: 2185-2189.
4. Aoyagi, S. and Archer, T.K. 2005. Modulating molecular chaperone HSP 90 functions through reversible acetylation. *Trends Cell Biol.* 15: 565-567.
5. Chen, B., Piel, W.H., Gui, L., Bruford, E. and Monteiro, A. 2005. The HSP 90 family of genes in the human genome: insights into their divergence and evolution. *Genomics* 86: 627-637.
6. Zhao, R. and Houry, W.A. 2005. HSP 90: a chaperone for protein folding and gene regulation. *Biochem. Cell Biol.* 83: 703-710.
7. Wegele, H., Wandinger, S.K., Schmid, A.B., Reinstein, J. and Buchner, J. 2005. Substrate transfer from the chaperone HSP 70 to HSP 90. *J. Mol. Biol.* 356: 802-811.
8. Allan, R.K., Mok, D., Ward, B.K. and Ratajczak, T. 2006. Modulation of chaperone function and cochaperone interaction by novobiocin in the C-terminal domain of HSP 90: evidence that coumarin antibiotics disrupt HSP 90 dimerization. *J. Biol. Chem.* 281: 7161-7171.

CHROMOSOMAL LOCATION

Genetic locus: Hsp90b1 (mouse) mapping to 10 C1.

PRODUCT

GRP 94 (m): 293T Lysate represents a lysate of mouse GRP 94 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

APPLICATIONS

GRP 94 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive GRP 94 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20 $^{\circ}$ C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.