



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# RGS11 (m): 293T Lysate: sc-123098

## BACKGROUND

Heterotrimeric G proteins function to relay information from cell surface receptors to various intracellular effectors. G proteins comprise  $\alpha$ ,  $\beta$  and  $\gamma$  subunits, and following activation the  $\alpha$  subunit binds GTP and dissociates from the  $\beta\gamma$  complex. A large group of proteins have been identified as GTPase-activating proteins (GAPs), including the RGS (regulator of G protein signaling) family, which serve to deactivate specific  $G_{\alpha}$  isoforms by increasing the rate at which they convert GTP to GDP. A subfamily of RGS proteins expressed in the central nervous system contain, in addition to the highly conserved RGS domain, a characteristic GGL domain, or G protein  $\gamma$  subunit-like domain, which mediates binding to  $G_{\beta 5}$  subunits. This subfamily, which includes RGS6, RGS7, RGS9 and RGS11, associates with  $G_{\beta 5}$  to form active GAP complexes that are predominantly localized to the cytosol. RGS/ $\beta 5$  complexes preferentially target  $G_{\alpha o}$  subunit for hydrolysis and inhibit  $G_{\beta 1\gamma 2}$ -mediated activation of phospholipase C.

## REFERENCES

1. Conklin, B.R. and Bourne, H.R. 1993. Structural elements of  $G_{\alpha}$  subunits that interact with  $G_{\beta\gamma}$  receptors, and effectors. *Cell* 73: 631-641.
2. Snow, B.E., Krumins, A.M., Brothers, G.M., Lee, S.F., Wall, M.A., Chung, S., Mangion, J., Arya, S., Gilman, A.G. and Siderovski, D.P. 1998. A G protein  $\gamma$  subunit-like domain shared between RGS11 and other RGS proteins specifies binding to  $G_{\beta 5}$  subunits. *Proc. Natl. Acad. Sci. USA* 95: 13307-13312.
3. Thomas, E.A., Danielson, P.E. and Sutcliffe, J.G. 1998. RGS9: a regulator of G-protein signalling with specific expression in rat and mouse striatum. *J. Neurosci. Res.* 52: 118-124.
4. Guan, K.L. and Han, M. 1999. A G-protein signaling network mediated by an RGS protein. *Genes Dev.* 13: 1763-1767.
5. Hepler, J.R. 1999. Emerging roles for RGS proteins in cell signalling. *Trends Pharmacol. Sci.* 20: 376-382.
6. Posner, B.A., Gilman, A.G. and Harris, B.A. 1999. Regulators of G protein signaling 6 and 7. Purification of complexes with  $G_{\beta 5}$  and assessment of their effects on G protein-mediated signalin pathways. *J. Biol. Chem.* 274: 31087-31093.
7. Snow, B.E., Betts, J.L., Mangion, J., Sondek, J. and Siderovski, D.P. 1999. Fidelity of G protein  $\beta$  subunit association by the G protein  $\gamma$  subunit-like domains of RGS6, RGS7, and RGS11. *Proc. Natl. Acad. Sci. USA* 96: 6489-6494.

## CHROMOSOMAL LOCATION

Genetic locus: Rgs11 (mouse) mapping to 17 A3.3.

## PRODUCT

RGS11 (m): 293T Lysate represents a lysate of mouse RGS11 transfected 293T cells and is provided as 100  $\mu$ g protein in 200  $\mu$ l SDS-PAGE buffer.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## APPLICATIONS

RGS11 (m): 293T Lysate is suitable as a Western Blotting positive control for mouse reactive RGS11 antibodies. Recommended use: 10-20  $\mu$ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

## STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.