



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# PDI (h): 293 Lysate: sc-158832

## BACKGROUND

Oxidoreductase-protein disulfide isomerase (PDI) is a homodimer consisting of subunits that catalyzes thiol-disulfide exchange, mediates folding of newly synthesized proteins and functions as a molecular chaperone. PDI localizes to the lumen of the endoplasmic reticulum (ER), where in conjunction with folding-helper proteins, such as immunoglobulin heavy chain binding protein (BiP), mediates tertiary and quaternary protein-processing. Cell surface PDI induces sulfhydryl-mediated conformational changes in integrin-mediated adhesion receptor-ligand interactions, thereby regulating integrin responses and cell adhesion. Additionally, PDI functions as a subunit of two more complex enzyme systems: the prolyl-4-hydroxylase and the triacylglycerol transfer proteins.

## REFERENCES

- Burgess, J.K., Hotchkiss, K.A., Suter, C., Dudman, N.P., Szollosi, J., Chesterman, C.N., Chong, B.H. and Hogg, P.J. 2000. Physical proximity and functional association of glycoprotein 1b $\alpha$  and protein-disulfide isomerase on the platelet plasma membrane. *J. Biol. Chem.* 275: 9758-9766.
- Klappa, P., Koivunen, P., Pirneskoski, A., Karvonen, P., Ruddock, L.W., Kivirikko, K.I. and Freedman, R.B. 2000. Mutations that destabilize the  $\alpha'$  domain of human protein-disulfide isomerase indirectly affect peptide binding. *J. Biol. Chem.* 275: 13213-13218.
- Mayer, M., Kies, U., Kammermeier, R. and Buchner, J. 2000. BiP and PDI cooperate in the oxidative folding of antibodies *in vitro*. *J. Biol. Chem.* 275: 29421-29425.
- Lahav, J., Gofer-Dadosh, N., Luboshitz, J., Hess, O. and Shaklai, M. 2000. Protein disulfide isomerase mediates Integrin-dependent adhesion. *FEBS Lett.* 475: 89-92.
- Maattanen, P., Kozlov, G., Gehring, K. and Thomas, D.Y. 2006. ERp57 and PDI: multifunctional protein disulfide isomerases with similar domain architectures but differing substrate-partner associations. *Biochem. Cell Biol.* 84: 881-889.
- Nuss, J.E., Choksi, K.B., DeFord, J.H. and Papaconstantinou, J. 2008. Decreased enzyme activities of chaperones PDI and BiP in aged mouse livers. *Biochem. Biophys. Res. Commun.* 365: 355-361.
- Appenzeller-Herzog, C. and Ellgaard, L. 2008. The human PDI family: versatility packed into a single fold. *Biochim. Biophys. Acta* 1783: 535-548.
- Shnyder, S.D., Mangum, J.E. and Hubbard, M.J. 2008. Triplex profiling of functionally distinct chaperones (ERp29/PDI/BiP) reveals marked heterogeneity of the endoplasmic reticulum proteome in cancer. *J. Proteome Res.* 7: 3364-3372.
- Hashimoto, S., Okada, K. and Imaoka, S. 2008. Interaction between bisphenol derivatives and protein disulfide isomerase (PDI) and inhibition of PDI functions: requirement of chemical structure for binding to PDI. *J. Biochem.* 144: 335-342.

## CHROMOSOMAL LOCATION

Genetic locus: P4HB (human) mapping to 17q25.3.

## PRODUCT

PDI (h): 293 Lysate represents a lysate of human PDI transfected 293 cells and is provided as 100  $\mu$ g protein in 200  $\mu$ l SDS-PAGE buffer.

## APPLICATIONS

PDI (h): 293 Lysate is suitable as a Western Blotting positive control for human reactive PDI antibodies. Recommended use: 10-20  $\mu$ l per lane.

Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

## STORAGE

Store at -20 $^{\circ}$  C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.