



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# DDC (h): 293T Lysate: sc-170280

## BACKGROUND

DOPA decarboxylase (DDC), also designated aromatic-L-amino-acid decarboxylase (AADC) belongs to the group II decarboxylase family of proteins. DDC, which can form a homodimer, is an important protein in the catecholamine biosynthesis pathway. DDC acts as a catalyst in the decarboxylation of L-5-hydroxytryptophan to Serotonin, L-3,4-dihydroxyphenylalanine (DOPA) to dopamine and L-tryptophan to tryptamine. Defects in the gene encoding for DDC may cause the autosomal recessive disorder AADC deficiency. AADC deficiency is an early onset inborn error in neurotransmitter metabolism which can lead to catecholamine and Serotonin deficiency. This causes poor feeding, psychomotor and developmental delays, lethargy, ptosis, gastrointestinal disturbances and hypothermia.

## REFERENCES

1. Scherer, L.J., McPherson, J.D., Wasmuth, J.J. and Marsh, J.L. 1992. Human DOPA decarboxylase: localization to human chromosome 7p11 and characterization of hepatic cDNAs. *Genomics* 13: 469-471.
2. Sumi-Ichinose, C., Ichinose, H., Takahashi, E., Hori, T. and Nagatsu, T. 1992. Molecular cloning of genomic DNA and chromosomal assignment of the gene for human aromatic L-amino acid decarboxylase, the enzyme for catecholamine and Serotonin biosynthesis. *Biochemistry* 31: 2229-2238.
3. Craig, S.P., Thai, A.L., Weber, M. and Craig, I.W. 1992. Localization of the gene for human aromatic L-amino acid decarboxylase (DDC) to chromosome 7p13 → p11 by *in situ* hybridisation. *Cytogenet. Cell Genet.* 61: 114-116.
4. Le Van Thai, A., Coste, E., Allen, J.M., Palmiter, R.D. and Weber, M.J. 1993. Identification of a neuron-specific promoter of human aromatic L-amino acid decarboxylase gene. *Brain Res. Mol. Brain Res.* 17: 227-238.
5. Vassilacopoulou, D., Sideris, D.C., Vassiliou, A.G. and Fragoulis, E.G. 2004. Identification and characterization of a novel form of the human L-DOPA decarboxylase mRNA. *Neurochem. Res.* 29: 1817-1823.
6. Chang, Y.T., Sharma, R., Marsh, J.L., McPherson, J.D., Bedell, J.A., Knust, A., Bräutigam, C., Hoffmann, G.F. and Hyland, K. 2004. Levodopa-responsive aromatic L-amino acid decarboxylase deficiency. *Ann. Neurol.* 55: 435-438.
7. Ma, J.Z., Beuten, J., Payne, T.J., Dupont, R.T., Elston, R.C. and Li, M.D. 2005. Haplotype analysis indicates an association between the DOPA decarboxylase (DDC) gene and nicotine dependence. *Hum. Mol. Genet.* 14: 1691-1698.

## CHROMOSOMAL LOCATION

Genetic locus: DDC (human) mapping to 7p12.1.

## PRODUCT

DDC (h): 293T Lysate represents a lysate of human DDC transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

## STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## APPLICATIONS

DDC (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive DDC antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.