

# Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

#### SANTA CRUZ BIOTECHNOLOGY, INC.

## caspase-4 (h2): 293T Lysate: sc-175856



#### BACKGROUND

A unique family of cysteine proteases has been described that differs in sequence, structure and substrate specificity from any previously described protease family. This family, termed CED-3/caspase-1, is comprised of caspase-1, caspase-2, caspase-3, caspase-4, caspase-6, caspase-7 (also designated Mch3, ICE-LAP3 or CMH-1), caspase-9 and caspase-10. CED-3/ caspase-1 family members function as key components of the apoptotic machinery and act to destroy specific target proteins which are critical to cellular longevity. Poly(ADP-ribose) polymerase plays an integral role in surveying for DNA mutations and double strand breaks. Caspase-3, caspase-7 and caspase-9, but not caspase-1, have been shown to cleave the nuclear protein PARP into an apoptotic fragment. Caspase-6, but not caspase-3, has been shown to cleave the nuclear lamins which are critical to maintaining the integrity of the nuclear envelope and cellular morphology. Caspase-10 has been shown to activate caspase-3 and caspase-7 in response to apoptotic stimuli.

#### REFERENCES

- 1. Lindahl, T., et al. 1995. Post-translational modification of poly (ADP-ribose) polymerase induced by DNA strand breaks. Trends Biochem. Sci. 20: 405-411.
- 2. Duan, H., et al. 1996. ICE-LAP3, a novel mammalian homologue of the Caenorhabditis elegans cell death protein CED-3 is activated during FASand tumor necrosis factor-induced apoptosis. J. Biol. Chem. 271: 1621-1625.
- 3. Duan, H., et al. 1996. ICE-LAP6, a novel member of the ICE/CED-3 gene family, is activated by the cytotoxic T cell protease granzyme B. J. Biol. Chem. 271: 16720-16724.
- 4. Fernandes-Alnemri, T.F., et al. 1996. In vitro activation of CPP32 and Mch3 by Mch4, a novel human apoptotic cysteine protease containing two FADD-like domains, Proc. Natl. Acad. Sci. USA 93: 7464-7469.
- 5. Simbulan-Rosenthal, C.M., et al. 1996. The expression of poly(ADP-ribose) polymerase during differentiation-linked DNA replication complex. Biochemistry 35: 11622-11633.
- 6. Casciola-Rosen, L., et al. 1996. Apopain/CPP32 cleaves proteins that are essential for cellular repair: a fundamental principle of apoptotic death. J. Exp. Med. 183: 1957-1964.
- 7. Takahashi, A., et al. 1996. Cleavage of lamin A by Mch2 a but not CPP32: multiple interleukin 1 β-converting enzyme-related proteases with distinct substrate recognition properties are active in apoptosis. Proc. Natl. Acad. Sci. USA 93: 8395-8400.
- 8. Kamada, S., et al. 1997. Involvement of caspase-4(-like) protease in Fasmediated apoptotic pathway. Oncogene 15: 285-290.
- 9. Hu, S., et al. 1998. Caspase-14 is a novel developmentally regulated protease. J. Biol. Chem. 273: 29648-29653.

#### **STORAGE**

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

#### CHROMOSOMAL LOCATION

Genetic locus: CASP4 (human) mapping to 11q22.3.

#### PRODUCT

caspase-4 (h2): 293T Lysate represents a lysate of human caspase-4 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

#### **APPLICATIONS**

caspase-4 (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive caspase-4 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

#### **RESEARCH USE**

For research use only, not for use in diagnostic procedures.

#### **PROTOCOLS**

See our web site at www.scbt.com for detailed protocols and support products.