



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-BCA-225, clone Cu-18

Clone no. Cu-18

MONOSAN Ready To Use

Product name	Mouse anti-BCA-225, clone Cu-18
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Polyclonal
Clone number	Cu-18
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-BCA-225, clone Cu-18

Clone no. Cu-18

MONOSAN Ready To Use

Additional info

The antibody recognizes a human breast carcinoma associated glycoprotein BCA-225 (220-225kD).

This protein differs in size and distribution from other breast carcinoma antigens. Anti-BCA-225

primary antibody labels breast cancer antigen 225 (BCA-225) in primary and metastatic breast

carcinoma. BCA-225 was first identified in T47D breast carcinoma cells, but its expression in other

carcinomas such as lung, kidney, ovary and endometrium has

References

1. Mesa-Tejada R, et al. Am J Pathol; 1988 130:305-14
2. Loy TS, et al. Am J Clin Pathol. 1991; 96:326-9
3. Brown RW, et al. Am J Clin Pathol. 1997; 107:12-9
4. Ansai S, et al. J Dermatol. 1994; 21:20-4
5. Sinard, JH. Arch Ophthalmol. 1999; 117:776-83

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES