



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Cytokeratin 8, clone 35betaH11

Clone no. 35betaH11

MONOSAN Ready To Use

Product name	Mouse anti-Cytokeratin 8, clone 35betaH11
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgM
Clonality	Monoclonal
Clone number	35betaH11
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cytokeratin 8, clone 35betaH11

Clone no. 35betaH11

MONOSAN Ready To Use

Additional info

Cytokeratin 8, a member of the Type II family of cytokeratins, is typically expressed in simple epithelium. The dimerization of cytokeratin 8 with cytokeratin 18 (labeled by 35betaH11) in the cytoplasm of simple epithelial cells allows for the formation of an intermediate filament cytoskeletal framework. This structure plays a role in the maintenance of cellular structural integrity and also functions in promoting signal transduction and cellular differentiation processes. Additionally, the presence of cytokeratin 8 has been detected in neoplastic epithelia, including glandular epithelium that can be found in prostate carcinoma. Positive immunoreactivity with anti-cytokeratin 8 is a useful indicator for the identification of normal and neoplastic epithelial tissues.

References

1. Battifora, H. Am J Surg Pathol 1988;12:24
2. Gown, AM, et al. Am J Clin Pathol 1985;84:413
3. Ljung G, et al. Prostate. 1997; 31:91-7
4. Murata T, et al. Pathol Res Pract. 1993; 189:888-93
5. Moll R, et al. Histochem Cell Biol. 2008; 129:705-33

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES