



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet

MON-RTU1095

MONOSAN[®]

Mouse anti-EMA, clone E29

Clone no. E29

MONOSAN Ready To Use

Product name	Mouse anti-EMA, clone E29
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a-k, IgG1-k
Clonality	Monoclonal
Clone number	E29
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-EMA, clone E29

Clone no. E29

MONOSAN Ready To Use

Additional info

The antibody is a useful marker for staining many carcinomas. It stains normal and neoplastic cells from various tissues, including mammary epithelium, sweat glands and colorectal carcinoma. Hepatocellular carcinoma, adrenal carcinoma and embryonal carcinomas are consistently EMA negative, so keratin positivity with negative EMA favors one of these tumors. EMA is frequently positive in meningioma, which can be useful when distinguishing it from other intracranial neoplasms such as schwannomas.

References

1. Dearnaly DP, et al. Lancet 1983; 1271-4
2. Attanoos RL, et al. Histopathology. 2003; 43:231-8
3. Beer TW, et al. Histopathology. 2000; 37:218-23
4. Lee JS, et al. Acta Cytol. 1996; 40:631-6
5. Fraga M, et al. Am J Clin Pathol. 1995; 103:82-9

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES