



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Epstein Barr Virus, clone CS1-4

Clone no. CS1-4

MONOSAN Ready To Use

---

Product name	Mouse anti-Epstein Barr Virus, clone CS1-4
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	CS1-4
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Epstein Barr Virus, clone CS1-4

Clone no. CS1-4

MONOSAN Ready To Use

---

**Additional info**

The antibody 4 targets the 60 kDa latent membrane protein (LMP-1) encoded by the BNLF1 gene of the Epstein-Barr virus. There is cross-reactivity with Reed Sternberg cells of Hodgkin's disease. The EpsteinBarr virus is an important as a cause of Infectious mononucleosis and has been associated with oral carcinomas.

**References**

1. Murray PG et al. J Pathol. 166: 1-5 (1992)
2. Jarrett RF et al. Blood 78:1-10 (1991)
3. - Paulesen G et al. Lancet. 337: 320-322 (1991)
4. Silverberg GS et al Principles and Practice of Surgical Pathology and Cytopatho
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**