



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Inhibin, alpha, clone R1

Clone no. R1

MONOSAN Ready To Use

Product name	Mouse anti-Inhibin, alpha, clone R1
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a
Clonality	Monoclonal
Clone number	R1
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Inhibin, alpha, clone R1

Clone no. R1

MONOSAN Ready To Use

Additional info

Inhibin is a peptide hormone that inhibits FSH secretion from the pituitary. Inhibin is a dimer that consists of an alpha and beta subunit. In normal tissue, anti-inhibin alpha labels granulosa cells of the ovary, sertoli and leydig cells of the testis, and the zona reticularis of the adrenal cortex. Anti-inhibin alpha has demonstrated utility in the identification of sex cord stromal tumors and adrenal cortical tumors.

References

1. Stewart CJ, et al. Histopathology. 1997; 31:67-74
2. McCluggage WG, et al. J Clin Pathol. 1998; 51: 114-6
3. Kommos F, et al. Mod Pathol. 1998; 11:656-64
4. Guerrieri C, et al. Int J Gynecol Pathol. 1998; 17:266-71
5. McCluggage WG. J Clin Pathol. 2000; 53:327-34

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES