



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Polyclonal antibody against LH

Clone no. -

MONOSAN Ready To Use

| | |
|---------------------------|---|
| Product name | Polyclonal antibody against LH |
| Host | Rabbit |
| Applications | IHC-P |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | - |
| Clonality | Monoclonal |
| Clone number | - |
| Size | 7 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against LH

Clone no. -

MONOSAN Ready To Use

Additional info

Luteinizing hormone (LH) is a heterodimeric glycoprotein produced by gonadotropic cells of the pituitary gland. Anti-LH is a useful marker to aid in the classification of pituitary tumors and the study of pituitary disease.

References

1. Sano T, et al. Virchows Arch A Pathol Anat Histopathol. 1990; 417:361-7
2. Felix I, et al. Hum Pathol. 1991; 22:719-21
3. Saccomanno K, et al. J Clin Endocrinol Metab. 1994; 78:1103-7
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES