



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-MART-1 + Tyrosinase, clone M2-7C10 & T311

Clone no. M2-7C10 &

MONOSAN Ready To Use

Product name	Mouse anti-MART-1 + Tyrosinase, clone M2-7C10 & T311
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2b/IgG2a
Clonality	Monoclonal
Clone number	M2-7C10 & T311
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MART-1 + Tyrosinase, clone M2-7C10 & T311

Clone no. M2-7C10 & T311

MONOSAN Ready To Use

Additional info

MART-1 (also known as Melan A) is a melanocyte differentiation antigen. It is present in melanocytes of normal skin, retina, nevi and in the majority of melanomas. Tyrosinase is an enzyme integral in the process of melanin synthesis, and found in most malignant melanomas. Therefore, this cocktail is useful for the identification of melanomas and melanocytic lesions

References

1. Orchard G. Br J Biomed Sci. 2002; 59:196-202
2. de Vries TJ, et al. J Pathol. 2001; 193:13-20
3. Shidham VB, et al. BMC Cancer. 2003; 3:15
4. Hoang MP, et al. J Cutan Pathol. 2001; 28:400-6
5. Ordóñez NG. Hum Pathol. 2014; 45:191-205

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES