



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Mouse anti-Melanoma, clone HMB-45

Clone no. HMB-45

MONOSAN Ready To Use

---

Product name	Mouse anti-Melanoma, clone HMB-45
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	HMB-45
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Melanoma, clone HMB-45

Clone no. HMB-45

MONOSAN Ready To Use

---

**Additional info**

Anti-HMB45 is a useful melanoma immunohistochemical marker that reacts with antigens present on immature melanosomes. Anti-HMB45 is useful for identifying amelanotic melanoma from other neoplastic lesions with similar morphology.

**References**

1. Gown AM, et al. Am J Pathol. 1986; 123:195-203
2. Wick MR, et al. Arch Pathol Lab Med. 1988; 112:616-20
3. Abrahamsen HN, et al. Cancer. 2004; 100:1683-91
4. Vaggelli L, et al. Tumori. 2000; 86:346-8
5. Baisden BL, et al. Am J Surg Pathol. 2000; 24:1140-6

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES