



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-WT1, clone 6F-H2

Clone no. 6F-H2

MONOSAN Ready To Use

Product name	Mouse anti-WT1, clone 6F-H2
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	6F-H2
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-WT1, clone 6F-H2

Clone no. 6F-H2

MONOSAN Ready To Use

Additional info

Wilms' tumor 1 protein (WT1) is a zinc finger transcription factor, normally expressed in tissues of mesodermal origin. The Wilms' tumor gene encodes a protein that functions as a tumor suppressor gene. WT1 is detected in tumor cells of Wilms' Tumor (also known as nephroblastoma) and mesothelioma. Additionally, WT1 expression has been found in ovarian serous carcinomas and some breast carcinomas.

References

1. Charles AK; Moore IE; Berry PJ. Histopathology ; 30(4):312-4 (1997)
2. Ordonez NG. Am J Surg Pathol 24(4):598-606, (2000)
3. Foster MR, et al. Arch Pathol Lab Med; 125:1316-20 (2001)
4. Nakatsuka S, et al. Mod Pathol; 19:804-14 (2006)
5. May RJ, et al. Clin Cancer Res.; 13: 4547-55 (2007)

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES