



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Polyclonal antibody against C4d

Clone no. -

MONOSAN Ready To Use

Product name	Polyclonal antibody against C4d
Host	Rabbit
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	-
Size	7 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against C4d

Clone no. -

MONOSAN Ready To Use

Additional info

C4d is a stable split product remnant of classical complement activation which becomes covalently bound to endothelium and basement membrane, after induction of the classical antibody-induced pathway. As an established marker of antibody-mediated acute renal allograft rejection and its proclivity for endothelium, this component can be detected in peritubular capillaries in both chronic renal allograft rejection as well as hyperacute rejection, acute vascular rejection, acute cellular rejection, and borderline rejection. It has been shown to be a significant predictor of transplant kidney graft survival and is an aid in treating acute rejection.

References

1. Jianghua C, et al. Clin Transplant. 2005; 19:785-91
2. Kayler LK, et al. Transplantation. 2008; 85:813-20
3. Ranjan P, et al. Nephrol Dial Transplant .2008; 23:1735-41
4. Seemayer CA, et al. Nephrol Dial Transplant. 2007; 22:568-76
5. Bouron-Dal Soglio D, et al. Hum Pathol. 2008; 39:1103-10

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES