



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44

MONOSAN Ready To Use

---

Product name	Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:15-1:30)
Species reactivity	human
Conjugate	-
Immunogen	DEAU cell line established from a large cell lymphoma
Isotype	IgMk
Clonality	Monoclonal
Clone number	DBA.44
Size	1 ml
Concentration	lot specific
Format	-
Storage buffer	Tissue culture supernatant with sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44

MONOSAN Ready To Use

---

**Additional info**

This antibody reacts with a subset of B-lymphocytes localized in the follicular mantle zone. It reacts with 97% of hairy cell leukemia cases. This antibody shows strong positive staining of about 35% of cases of high grade B cell lymphomas. Positive control Tonsil.

**References**

1. Al Saati et al. 1989;74:2476
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES