



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Rabbit anti-PROGESTERON RECEPTOR, clone ZR4 (Monoclonal)

Clone no. ZR4

MONOSAN

---

Product name	Rabbit anti-PROGESTERON RECEPTOR, clone ZR4 (Monoclonal)
Host	Rabbit
Applications	IHC-P (1:200-1:400)
Species reactivity	Human
Conjugate	-
Immunogen	Recombinant protein encoding human progesterone receptor 412-526 aa
Isotype	IgG
Clonality	Monoclonal
Clone number	ZR4
Size	1 ml
Concentration	n/a
Format	Concentrate
Storage buffer	Tissue culture supernatant with 0.2% BSA and 15mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rabbit anti-PROGESTERON RECEPTOR, clone ZR4 (Monoclonal)

Clone no. ZR4

MONOSAN

**Additional info**

The progesterone receptor (PgR) is an estrogen regulated protein. It has been proposed that expression of PgR determination indicates a responsive estrogen receptor (ER) pathway, and therefore, may predict likely response to endocrine therapy in human breast cancer. A number of studies have shown that PgR determination provides supplementary information to ER, both in predicting response to endocrine therapy and estimating survival. PgR has proved superior to ER as a prognostic indicator in some studies.

Pretreatment: Heat induced epitope retrieval in 10 mM citrate buffer, pH6.0 for 20 minutes is required for IHC staining on formalin-fixed, paraffin embedded tissue sections. Note: Dilution of the antibody in 10% normal goat serum followed by a goat anti-rabbit secondary antibody-based detection is recommended. Control tissue Breast, breast ductal carcinoma. Staining Nuclear

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**