



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet MON10223

MONOSAN[®]

Mouse anti-E-CADHERIN, clone MX020 (Monoclonal)

Clone no. MX020

MONOSAN

Product name	Mouse anti-E-CADHERIN, clone MX020 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:100-1:200)
Species reactivity	Human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG
Clonality	Monoclonal
Clone number	MX020
Size	1 ml
Concentration	n/a
Format	Concentrate
Storage buffer	Tissue culture supernatant containing 15mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-E-CADHERIN, clone MX020 (Monoclonal)

Clone no. MX020

MONOSAN

Additional info

ECAD (syn, CD325), 120 kDA, chromosome 16q22.1 (CDH1 gene), is a critical regulator of epithelial junction formation. It interacts with the cytoskeleton through several associated proteins. The ECAD internal domain binds with alpha, beta, gamma and p120 catenin's to anchor the ECAD complex to the actin cytoskeleton of the cell. ECAD is expressed in virtually all epithelial cells except for adrenocortical cells. The expression in liver cells is weaker than in most other epithelia. ECAD is also expressed in melanocytes (adhering to squamous epithelial cells).

Pretreatment: Heat induced epitope retrieval in 10 mM citrate buffer , pH6.0 for 20 minutes is required for IHC staining on formalin-fixed, paraffin embedded tissue sections. Note: Dilution of the antibody in 10% normal goat serum followed by a goat anti-mouse secondary antibody-based detection is recommended. Control tissue Pancreas, lung adenocarcinoma, breast.

Staining Membranous

References

1. Ceyran A B, Şenol S, Şimşek B Ç, et al. Int journal of clinical and exp.pathol,201
2. Liu D, Huang C, Kameyama K, et al. The Annals of thoracic surgery, 2001, 71(3):
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES